
PATENTED TECHNOLOGY

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
ENGLISH KNEE BRACE EPI THELIUMFLEX® 01

INDICATIONS: Wearing the EPI THELIUMFLEX® 01 knee brace helps relieve patellar conditions (patellar syndrome, patellar instability, patellar tendonitis, etc.). In the event of a knee injury with residual effects (ligament weakness, osteoarthritis, etc.), it helps stabilise the joint when resuming sports activities.

FUNCTIONAL: Its patented EPI THELIUMFLEX® silicone patellar support ring combined with REFLEX ultra-effective compression fabric improves kneecap support and knee stability.

SPECIFIC: The knee brace works dynamically, without impeding movement. It can therefore be used preventively since it maintains muscle strength.

ERGONOMIC: Concentrated technology weighing less than 60g, the EPI THELIUMFLEX® 01 proprioceptive knee brace has been specially designed for sports. The EPI THELIUMFLEX® 01 knee brace is designed to fit either the right or the left knee. It does not slip and does not cause any discomfort at the back of the knee (popliteal fossa). Extra-thin, it can be worn under all types of clothing.

IMPORTANT! If you use calf compression sleeves (such as Booster), place the lower elastic band of the knee brace over the sleeve.

HOLD:
 - Do not apply moisturising products before using the knee brace to ensure optimum hold.
 - If you encounter any hold problems in extreme conditions (sports with sudden changes of speed, running on steep slopes or over long distances, high temperatures), we recommend that you wear a thigh sleeve over the knee brace.

IMPORTANT !
 - Place the knee brace on the thigh without stretching.
 - To guarantee good hold, the product must be washed after each use. At the very least, wash the silicone elastic bands with soap and water (without rubbing) to remove any dirt.



WASHING INSTRUCTIONS: Machine wash at 40°C or hand wash (without rubbing), leave to dry naturally. Do not iron.

CONTRAINDICATIONS: To ensure blood circulation is not impeded, do not use in the event of arteritis, severe varicose conditions or residual effects of phlebitis or oedema (lymphatic/venous) of the lower leg and/or thigh. Do not use with compression hosiery. Do not use on an open wound.

IMPORTANT! To avoid any risks of irritation generated by the silicone elastic bands, the knee brace should not be worn in static positions, particularly with the wearer lying down (sleep) or sitting (including bus or car journeys, etc.). This knee brace should be put on immediately before activity, to be worn during movement, during sports activities.

PRECAUTIONS FOR USE:
 • The product has undergone all the regulatory tests required by ISO standard 10993-1 relative to products that come into contact with the skin. However, if an incident occurs, discontinue use of the product immediately and inform MILLET Innovation.
 • If you experience pain, oedema, discomfort or any other abnormal reaction, stop using the product and consult a doctor.

COMPOSITION: 30% silicone, 70% polyamide/elastane fabric.
 MILLET Innovation cannot be held liable in cases where the manufacturer's recommendations have not been followed during use.

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GB HOW TO FIT THE KNEE BRACE

A) Unfold the elastic bands outwards.
 B) Fold the knee brace back on itself to slip it on.
 C) Position the kneecap in the centre of the "support ring".
 D) Unfold the upper elastic band, placing the knee brace on the thigh without tension.
 E) F) Unfold the lower elastic band and adjust your knee brace in the standing position.

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